

Special Soft Resin Denture Lining

About This Unique Material.

You have been prescribed a soft lining for the tissue surface of your denture, the portion next to the roof or floor of your mouth. The soft lining is a special compounded dental resin of unusual density and purity. At mouth temperature, this lining has the same resiliency and softness as normal mouth tissue, creating a cushion which dramatically improves the fit and comfort of your denture. This resin is tasteless, odourless, and stain resistant. With proper care and regular examinations by your dentist, it will give you better service than any other material currently available.

About Care & Cleaning.

You may clean your denture with most denture cleansers or toothpaste's, following the directions of the manufacturer and the instructions of your dentist. Regular cleaning is important for good oral hygiene. Before cleaning the denture, part-fill the wash basin with cold tap water and immerse the denture for a few minutes. Cold water will temporarily harden the soft lining; warm water will soften it. Using a soft toothbrush and cleanser of your choice, scrub the denture gently. Hold the denture above a basin part-filled with water. If the denture is accidentally dropped, the water will cushion the fall, preventing damage to the teeth. When washing a lower denture don't grasp it at both ends since too much pressure might cause it to break in the centre. Rinse the denture thoroughly in warm water. The soft lining will return to its original resiliency when the denture returns to mouth temperature.

Do not soak your denture in solutions containing bleach - this may cause the soft lining to harden and change colour!

